
Continuous Glucose Monitoring assessing the Clinical Impact of an Algorithm Driven Basal-Bolus Insulin Regimen in non-critically ill inpatients with T2DM

K. Donsa, K. M. Neubauer, J. K. Mader, B. Höll, S. Spat, P. Beck, F. Aberer, T. Augustin, J. Plank, T. R. Pieber, L. Schaupp

Medical University of Graz, Austria
Division of Endocrinology and Metabolism

JOANNEUM RESEARCH
Forschungsgesellschaft mbH

HEALTH – Institute for
Biomedicine and Health Sciences

Klaus Donsa
21.02.2015



Disclosure

X	No, nothing to disclose
	Yes, please specify:

<i>Company Name</i>	<i>Honoraria/ Expenses</i>	<i>Consulting/ Advisory Board</i>	<i>Funded Research</i>	<i>Royalties/ Patent</i>	<i>Stock Options</i>	<i>Ownership/ Equity Position</i>	<i>Employee</i>	<i>Other (please specify)</i>

Background

Inpatient glycemic management

3

- Improvement in glycemic control results in lower rates of hospital complications (Umpierrez et al. 2012, Murad et al. 2012)
- Guideline recommendations (ADA 2014)
 - Standardized insulin order set
 - Scheduled insulin
 - Pre-meal BG values <140 mg/dl
 - Four capillary BG measurements
- Four point BG profile provides only a sequence of snapshots (Joseph et al. 2009, Thabit 2014)

Background

Inpatient glycemic management

4


- CGM is used
 - ...to detect patterns of responsiveness to therapeutic efforts (Gomez 2014)
 - ...as a tool to assess the outcome of clinical trials (Beck 2012)

Objective


5

- To
 - ...assess the clinical impact of an algorithm driven basal-bolus-insulin regimen
 - ...compare CGM to routine blood glucose measurements in the hospital

Trial setting



- BG POCT (Roche, Accu Chek)
- CGM (Medtronic, iPro2)
- Insulin injection
- Bolus insulin (meal)
- Bolus insulin (supplement)
- Basal insulin



Patient population

- **General ward.** Division for Endocrinology and Metabolism (Department of Internal Medicine, Medical University Graz, Austria)
- **Main inclusion criteria**
 - T2DM
 - Age between 18 and 90 years
 - Any antihyperglycemic therapy
- **Main exclusion criteria**
 - T1DM
 - Pregnancy
 - Terminal illness

Characteristics	
n	84
Gender, f/m	32/52
Age (years)	68.5 ± 10.3
BMI (kg/m ²)	31.0 ± 6.6
Serum creatinine (mg/dl)	1.5 ± 1.0
HbA1c (mmol/mol)	72.8 ± 28.2
Diabetes duration (years)	15.1 ± 11.1
Admission diagnosis (%)	
Cardiovascular diseases	31.3
Endocrine disorder	25.0
Infectious diseases	30.2
Gastrointestinal disease	1.0
Other	12.5

Comparison between CGM & BG

8

Glucose Profile	Overall		First Day		Last Day		Test
	CGM	BG	CGM	BG	CGM	BG	
Glucose values (n)	140,424	2,066	23,686	351	23,301	323	
Mean glucose (mg/dl)	145.9	149.6	156.8	160.4	142.3	145.9	<0.02
SD (mg/dl)	55.9	57.7	63.1	64.9	52.3	54.1	
Coefficient of variation CV (%)	37.7	38.6	39.6	40.6	36.9	36.8	<0.03/ 0.05
<50 mg/dl (%)	0.3	0.2	0.5	0.0	0.4	0.0	
<70 mg/dl (%)	2.4	2.7	2.6	1.7	2.8	1.2	>0.2
70-180 mg/dl (%)	75.8	73.2	67.7	67.2	77.5	78.6	<0.04
>180 mg/dl (%)	21.9	24.1	29.7	31.1	19.7	20.1	<0.04
>350 mg/dl (%)	0.5	0.6	0.9	1.1	0.4	0.6	

Comparison between CGM & BG

9

Glucose Profile	Overall		First Day		Last Day		Test
	CGM	BG	CGM	BG	CGM	BG	p-Value
Glucose values (n)	140,424	2,066	23,686	351	23,301	323	
Mean glucose (mg/dl)	145.9	149.6	156.8	160.4	142.3	145.9	<0.02
SD (mg/dl)	55.9	57.7	63.1	64.9	52.3	54.1	
Coefficient of variation CV (%)	37.7	38.6	39.6	40.6	36.9	36.8	<0.03/ 0.05
<50 mg/dl (%)	0.3	0.2	0.5	0.0	0.4	0.0	
<70 mg/dl (%)	2.4	2.7	2.6	1.7	2.8	1.2	
70-180 mg/dl (%)	75.8	73.2	67.7	67.2	77.5	78.6	<0.04
>180 mg/dl (%)	5.6	5.9	8.6	10.0	3.7	3.4	
>350 mg/dl (%)	0.5	0.6	0.9	1.1	0.4	0.6	

Significant reduction of the CV

Comparison between CGM & BG

10


Glucose Profile	Overall		First Day		Last Day		Test
	CGM	BG	CGM	BG	CGM	BG	p-Value
Glucose values (n)	140,424	2,066	23,686	351	23,301	323	
Mean glucose (mg/dl)	145.9	149.6	156.8	160.4	142.3	145.9	<0.02
SD (mg/dl)	55.9	57.7	63.1	64.9	52.3	54.1	
Coefficient of variation CV (%)	37.7	38.6	39.6	40.6	36.9	36.8	<0.03/ 0.05
<50 mg/dl (%)	0.3	0.2	0.5	0.0	0.4	0.0	
<70 mg/dl (%)	2.4	2.7	2.6	1.7	2.8	1.2	>0.2
70-180 mg/dl (%)	75.8	73.2	67.7	67.2	77.5	78.6	<0.04
>180 mg/dl (%)	21.9	24.1	29.7	31.1	19.7	20.1	<0.04
>350 mg/dl (%)	0.5	0.6	0.9	1.1	0.4	0.6	

Significant increase in the range 70-180 mg/dl

Glucose profile

First full day with CGM


11



Glucose profile

Last full day with CGM


12



Clinical impact

First vs. last full day with CGM

13



Conclusion

14

- High consistency between CGM and BG measurements
- CGM describes the overall daily routine in more detail
- Assessing the Basal-Bolus Therapy:
 - Lower glucose variations at the end of the hospital stay
 - Therapy yielded a higher percentage of BG values in the range 70-180 mg/dl
- A trend towards detecting more low-glycemic events with CGM



Medical University of Graz



Thank you for your attention

DI Klaus Donsa, BSc

mail: klaus.donsa@joanneum.at

The study is supported by the European Commission
Project REACTION: FP7-248590



Medical University of Graz
Division of Endocrinology and
Metabolism
Auenbruggerplatz 15, 8036 Graz
www.medunigraz.at

JOANNEUM RESEARCH
Forschungsgesellschaft mbH
HEALTH – Institute for Biomedicine
and Health Sciences
Neue Stiftingtalstraße 2, 8010 Graz
www.joanneum.at/health

Dosing algorithm

A

initial TDD	TDD: 48 IU	yesterday lunch (IU)	yesterday dinner (IU)	yesterday bedtime (IU)	today breakfast (IU)
	basal insulin (TDD/2)	24			
	bolus insulin (TDD/2/3)	8	8		8
	BG (mmol/l)	9.2	7.7	7.9	10.6



		breakfast: 10.6 mmol/l			
dinner: 7.7 (mmol/l)	BG range (mmol/l)	3.9 – <5.6	5.6 - 7.8	>7.8 – 10.0	>10.0
	3.9 – <5.6	-10 %	± 0 %	± 0 %	± 0 %
	5.6 – 7.8	± 0 %	± 0 %	± 0 %	+ 10 %
	>7.8 – 10.0	± 0 %	± 0 %	+ 10 %	+ 10 %
	>10.0	± 0 %	+ 10 %	+ 10 %	+ 20 %
	if BG ≤ 3.9 mmol/l reduce next 24 h TDD by 20%				



new TDD	next 24 h TDD: 52 IU	today lunch (IU)	today dinner (IU)	today bedtime (IU)	tomorrow breakfast (IU)
	basal insulin (TDD/2)	26			
	bolus insulin (TDD/2/3)	9 (8+1)	8		9 (8+1)

Dosing algorithm

insulin sensitivity mmol/l	O sensitive	O normal	O resistant
< 3.9	no bolus insulin		
3.9 – <5.6	reduce suggested bolus insulin dose by 50%		
5.6 – 7.8	0	0	0
>7.8 – 10.0	2	4	6
>10.0 – 12.2	4	6	8
>12.2 – 14.4	6	8	10
>14.4 – 16.7	8	10	12
>16.7 – 19.4	10	12	14
>19.4 – 22.2	12	14	16
> 22.2	14	16	18